## CAD/CAM ABUTMENTS

<b>EN</b>	<b>Instruction for Use</b> <b>AoN Implants US</b> prosthetic components: CAD/CAM Abutments	<b>2 - 5</b>
<b>ES</b>	<b>Instrucciones de Uso</b> <b>AoN Implants US</b> componentes protésicos: Pilares CAD/CAM	<b>6 - 9</b>

- All AoN Implants CAD/CAM abutments are only intended to be used as straight abutments (i.e., 0° angulation).
- Compatible Implant Systems: K192015, AoN Implants Dental Implant System

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**Caution: Federal Law restricts this device to sale by or on the order of a licensed physician or health-care provider.**

**GENERAL INFORMATION:**

AoN Implants CAD/CAM Abutments are CAD-CAM abutments made from titanium bases (T-Base) with zirconia copings. Some of the T-Bases are CEREC compatible and some are Multi-Unit T-Bases. The titanium bases connect to the currently cleared AoN Implants dental implants. The connection platforms for each of these bases are compatible with our internal hex, external hex, and Conical implants.

Dentists can scan for the required dimensions and type of base and zirconia coping; these scans are sent to AoN Implants and AoN Implants will use CAD to design an appropriate abutment [T-Base and zirconia coping] to be used with the appropriate implant. AoN Implants will then mill the abutments using CAM at their GMP facility which is registered with FDA under FURLS. AoN Implants will also do the final assembly.

The following summary provides a snapshot of the devices.

- There are 3 different T-bases which are available in 3 different platforms. Each T-base is available in each platform.

Titanium Base	Platform Available		
	Internal Hex Platform	External Hex Platform	Standard Conical Platform
T-Base for CEREC	Yes	Yes	Yes
T-Base	Yes	Yes	Yes
MU T-Base	Yes	Yes	Yes


- The Zirconia Coping is applicable to all T-bases including MU and CEREC.

Zirconia Coping	T-Base Available		
	T-Base for CEREC	T-Base	MU T-Base
	Yes	Yes	Yes

- None of these T-Bases or Zirconia Copings have angulation, all angulations are 0 degrees.
- All these T-bases, MU T-Bases, and CEREC T-bases will be assembled to the zirconia coping by AoN Implants technicians in a GMP facility.
- These Straight Multi-Unit (MU) Abutments are not the previously cleared Multi-Unit (MU) Abutments, they are Titanium Bases (T-Base) intended to be used with a Zirconia Coping in the exact same way as the other T-Base/Zirconia Coping devices in this family. The MU T-Base does not contain angulation.
- The "zirconia" ceramic material that comprises the coping of the two-piece T-Base abutments are FDA cleared under K142987, DD Bio Z Blanks.
- AoN Implants will not provide an unfinished T-Base (without ceramic coping) to the end user for assembly. The T-Base CEREC abutments will only be manufactured and assembled by AoN Implants to ensure the proposed design parameters are being followed.
- The AoN Implants Dental Implant/Abutment System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the AoN Implants Dental Implant/Abutment System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Illustration	Description	Materials
	Dental Cement	K982394 (Prime & Bond) K040906 (Calibra)
	Zirconia Coping	K142987 (DD Bio Z)
	Titanium Base	K192015 (TI 6AL4V)
	Sterilization	K192015 (Autoclave)
	Computer Aided Design (CAD)	AoN Implants -validated
	Computer Aided Manufacturing (CAM) Milling	AoN Implants -validated
	Assembly	AoN Implants -validated

The crown copings validated for these CAD/CAM abutments should be within the following ranges: 7mm-13mm(W) x 7mm-13mm(D) x 11mm-19mm(H) with 0.5mm-5mm (wall thickness at base) coping.

(Premilled Blank) Customizable Abutment	
Description	The Customizable Abutment is an abutment with an internal hex or external hex platform that connects to the reference K192015 implants for use in edentulous situations
Restoration Technique	to be cemented on the coping
Characteristics	The abutment is intended for single use. It is fixated using a cement to the coping and securement screws to the implant. The prosthetic is not used with another prosthetic.
Family Type	T-base (CAD/CAM)
Implant Attachment	Direct Screwed
Prosthetic Attachment	K982394 Prime & Bond, K040906 Calibra
Material	Titanium
Sterilization	Autoclave
Angulation (°)	0
Anti-Rotation	Hex
Diameter	4 - 12 mm
Wall Thickness	0.5-5mm
Length	12mm
Post Height	4 - 7mm
Total Height	7-12mm
Gingival Height	0-3mm
Zirconia Coping Width	7mm-13mm
Zirconia Coping Depth	7mm-13mm
Zirconia Coping Height	11mm-19mm
Zirconia Coping Wall Thickness at Base	0.5mm-5mm

T-Base "S" (CEREC)	
Description	The T-base is an abutment with an internal hex, conical, or external hex platform that connects to the reference K192015 implants for use in edentulous situations
Restoration Technique	to be screwed on the implant
Characteristics	The abutment is intended for single use. It is fixated using a hex driver that turns and tightens the securement screws to the implant. The abutment is fitted with a zirconia coping.
Family Type	T-base (CAD/CAM)
Implant Attachment	Securement Screw
Prosthetic Attachment	K982394 Prime & Bond, K040906 Calibra
Material	ASTM F-136
Sterilization	Autoclave
Angulation (°)	0
Diameter	4.3 - 5 mm
Length	5.10 - 13mm
Shoulder Height	0.5 - 4.5mm
Connection Platform	External hex    Internal Hex    Conical
Platform Height	3.4 - 4.10mm    3.5mm    2.8mm
Gingival Height	0.5 - 1.5 - 3.0mm    0.5 - 1.5 - 3.0mm    0.7 - 1.5 - 3.0 - 4.5mm
Post Height	4.6mm    4.6mm    4.6mm
Internal Hex Width	2.5 - 2.7mm    2.43mm    2.20mm
Internal Hex Depth	0.7mm    1.7mm    3.5mm
Platform Taper	0 deg    45 deg    3 deg
Zirconia Coping Width	7mm-13mm
Zirconia Coping Depth	7mm-13mm
Zirconia Coping Height	11mm-19mm
Zirconia Coping Wall Thickness at Base	0.5mm-5mm

T-Base 3Pod	
Description	The T-Base 3Pod is an abutment with an internal hex, conical, or external hex platform that connects to the reference K192015 implants for use in edentulous situations
Restoration Technique	to be screwed on the implant
Characteristics	The abutment is intended for single use. It is fixated using a hex driver that turns and tightens the securement screws to the implant. The abutment is fitted with a zirconia coping.
Family Type	T-base (CAD/CAM)
Implant Attachment	Securement Screw
Prosthetic Attachment	K982394 Prime & Bond, K040906 Calibra
Material	ASTM F-136
Sterilization	Autoclave
Angulation (°)	0
Diameter	4.7mm
Length	4.5 - 12.4mm
Shoulder Height	0.5 - 4.5mm

T-Base 3Pod			
Connection Platform	External hex	Internal Hex	Conical
Platform Height	3.4 – 4.10mm	3.5mm	2.8mm
Gingival Height	0.5 – 1.5 – 3.0mm	0.5 – 1.5 – 3.0mm	0.7 – 1.5 – 3.0 – 4.5mm
Post Height	4.0mm	4.0mm	4.0mm
Internal Hex Width	2.5 – 2.7mm	2.43mm	2.20mm
Internal Hex Depth	0.7mm	1.7mm	3.5mm
Platform Taper	0 deg	45 deg	3 deg
Zirconia Coping Width	7mm-13mm		
Zirconia Coping Depth	7mm-13mm		
Zirconia Coping Height	11mm-19mm		
Zirconia Coping Wall Thickness at Base	0.5mm-5mm		

T-Base Abutment for Transepithelial Connector			
Description	The T-Base is an Abutment with the intention of use in conjunction with the Transepithelial Connector that connects to the reference K192015 implants for use in edentulous situations.		
Restoration Technique	to be screwed on transepithelial connector		
Characteristics	The abutment is intended for single use. It is fixated using a hex driver that turns and tightens the securement screws to the transepithelial connector. The abutment is fitted with a zirconia coping.		
Family Type	Transepithelial Connector (CAD/CAM)		
Implant Attachment	Securement Screw		
Prosthetic Attachment	K982394 Prime & Bond, K040906 Calibra		
Material	ASTM F-136		
Sterilization	Autoclave		
Angulation (°)	0		
Diameter	4.1 – 5.0mm		
Length	4.3 – 4.9mm		
Shoulder Height	0.30mm		
Connection Platform	Internal Hex	Conical	
Platform Height	3.4 – 5.0mm	3.4 – 5.0mm	
Zirconia Coping Width	7mm-13mm		
Zirconia Coping Depth	7mm-13mm		
Zirconia Coping Height	11mm-19mm		
Zirconia Coping Wall Thickness at Base	0.5mm-5mm		

T-Base 3Pod for Transepithelial Connector			
Description	The T-Base 3Pod is an Abutment with the intention of use in conjunction with the Transepithelial Connector that connects to the reference K192015 implants for use in edentulous situations.		
Restoration Technique	to be screwed on transepithelial connector		

T-Base 3Pod for Transepithelial Connector			
Characteristics	The abutment is intended for single use. It is fixated using a hex driver that turns and tightens the securement screws to the transepithelial connector. The abutment is fitted with a zirconia coping.		
Family Type	T-base (CAD/CAM)		
Implant Attachment	Securement Screw		
Prosthetic Attachment	K982394 Prime & Bond, K040906 Calibra		
Material	ASTM F-136		
Sterilization	Autoclave		
Angulation (°)	0		
Diameter	4.7 – 5.0mm		
Length	4.3 – 4.9mm		
Shoulder Height	0.30mm		
Connection Platform	Internal Hex	Conical	
Platform Height	3.4 – 5.0mm	3.4 – 5.0mm	
Zirconia Coping Width	7mm-13mm		
Zirconia Coping Depth	7mm-13mm		
Zirconia Coping Height	11mm-19mm		
Zirconia Coping Wall Thickness at Base	0.5mm-5mm		

T-Base 3Pod Multi Utility			
Description	The T-Base 3Pod Multi Utility is an Abutment with the intention of use in conjunction with the Multi Utility abutment that connects to the reference K192015 implants for use in edentulous situations.		
Restoration Technique	to be screwed on the Multi Utility		
Characteristics	The abutment is intended for single use. It is fixated using a hex driver that turns and tightens the securement screws to the Multi Utility. The abutment is fitted with a zirconia coping.		
Family Type	T-base (CAD/CAM)		
Implant Attachment	Securement Screw		
Prosthetic Attachment	K982394 Prime & Bond, K040906 Calibra		
Material	ASTM F-136		
Sterilization	Autoclave		
Angulation (°)	0		
Diameter	5mm		
Length	4.3mm		
Shoulder Height	0.3mm		
Connection Platform	External hex	Internal Hex	Conical
Platform Height	5mm	5mm	5mm
Gingival Height	0.5 – 1.5 – 3.0mm	0.5 – 1.5 – 3.0mm	0.7 – 1.5 – 3.0 – 4.5mm
Post Height	5.1 – 6.1 – 7.6mm	5.1 – 6.1 – 7.6mm	5.1 – 6.1 – 7.6mm
Internal Hex Width	2.5 – 2.7mm	2.43mm	2.20mm
Internal Hex Depth	0.7mm	1.7mm	3.5mm
Platform Taper	0 deg	45 deg	3 deg
Platform Height	5mm	5mm	5mm

T-Base 3Pod Multi Utility	
Zirconia Coping Width	7mm-13mm
Zirconia Coping Depth	7mm-13mm
Zirconia Coping Height	11mm-19mm
Zirconia Coping Wall Thickness at Base	0.5mm-5mm

### Indications for Use:

AoN Implants CAD/CAM Abutments are dental abutments placed onto a dental implant to provide support for dental prosthetic restorations. The abutments include: 1) Titanium abutment blanks with a pre-machined implant connection where the upper portion may be custom-milled in accordance with a patient-specific design using CAD/CAM techniques; and 2) Titanium bases with a pre-machined implant connection upon which a CAD/CAM designed superstructure may be fitted to complete a two-piece dental abutment. The abutments include an abutment screw for fixation to the underlying implant. The abutments may be used for single-unit (single-tooth) or multiple-unit (bridges and bars) restorations and are compatible for use with AoN Implants Dental Implants with internal hex-connection mating platform and conical connection platform diameters.

All digitally designed abutments and/or copings for use with AoN Implants CAD/CAM Abutments are intended to be sent to a AoN Implants-validated milling center for manufacture. AoN Implants abutments designed using CAD/CAM techniques must fulfill the AoN Implants allowable range of design parameters.

### Directions for Use

1. **Scan Body** must be selected according to compatible Implant system and connection platform.

#### Extra-oral Scanning process

Soft tissue and gingival mask should be removed from the model cast.

Screw the Scan Body (hand tightening), so that it perfectly seated on the implant or Multi Unit analog.

#### Intra-oral Scanning process

Before use, the Scan Body should be sterilized. Screw the Scan Body (hand tightening), so that it perfectly seated on the implant or Multi Unit.

For most of the scanners, there is no need to use spray, otherwise, apply a thin coat of scan spray from a 20 cm distance.

During the scanning of a bridge, all Scan Bodies should be scanned at the same time and not separately to eliminate inaccuracies that may occur during screwing and unscrewing.

The Scan Body and fixing screw are sterilizable and are intended for single use.

2. **T-Base** must be selected to be compatible with the original implant system and connection platform.

There are two types of T-Base:

- With anti-rotational lock (engaged) – intended for single tooth restoration.
- Without anti rotational lock (non-engaged) – intended for multiple teeth restoration.

The diameter of the T-Base must not be reduced e.g., by grinding. Shortening of the T-Base is not recommended.

The contact surfaces of the T-Base with the implant should not be sand-blasted or modified. Before sandblasting the upper part, the Ti-Base should be fastened in a laboratory analog to protect the implant interface area. Only the surfaces of the T-Base intended for gluing with a mesostructures needs to be sandblasted (50-120µm aluminum oxide, max. 2.0 bar) and then cleaned (with alcohol or steam).

Use resin cement as an adhesive to connect the T-Base and the prosthetic structure. (According to the glue manufacturer's instructions).

For easier handling during the gluing process, it is recommended to screw the T-Base into a lab implant or a polishing tool.

Cover the hex head of the abutment screw with wax.

Use resin cement as an adhesive to connect the T-Base and the prosthetic structure. (According to the glue manufacturer's instructions).

Push the prosthetic structure in as far as it will go. Make sure it latches into the rotation and position stops.

Remove excess glue immediately.

Remove residue with a rubber polisher after hardening.

The dentist must sterilize the T-Base and the mesostructured before inserting it in the patient's mouth.

Recommended final restoration tightening torques:

- 25 Ncm when fixing the bases on Internal Hex platform.
  - 25 Ncm when fixing the bases on External Hex platform.
  - 25 Ncm when fixing the bases on Conical platform.
3. **Scan Post** must be selected according to compatible implant system and connection platform.

#### Extra-oral Scanning process

Soft tissue and gingival mask should be removed from the model cast.

Screw the Scan Post (hand tightening), so that it is perfectly seated on the implant analog.

#### Intra-oral Scanning process

Before use, Scan Post should be sterilized. Screw the Scan Post (hand tightening), so that it is perfectly seated on the implant. Scan body should be mounted on the scan post for scanning.

During the scanning of a bridge, all Scan Posts should be scanned at the same time, and not separately to eliminate inaccuracies that may occur during screwing and unscrewing.

Scan Post must not be used for the final implant restoration.

Subject AoN Implants			
Abutment	CYROTH Compatible	IS FOUR Compatible	KA-ROTH Compatible
Transepithelial Connector	Yes	Yes	Yes
Overdenture Prosthetic	Yes	Yes	Yes

#### CONTRAINDICATIONS:

It is contraindicated placing abutments in the following patients:

1. Medically unfit for oral surgery procedure (on corticosteroids, or anticoagulants or bisphosphonates, and those receiving radiation of other immunosuppressive therapy).
2. Lactating or pregnant women are not candidates, nor are patients with abnormal laboratory values for BUN, creatinine, or serum calcium.
3. Patients with uncontrolled diabetes, cardiovascular disease, and Hypertension above 170/110 mm Hg.
4. Osteoporotic crush fractures, respiratory disease, thyroid or parathyroid as well as patients with diagnosed malignancy disease or unexplained lump or masses in the head or neck.
5. Patients with uncontrolled disease such as Hemophilia, Granulocytopenia or other bleeding problems, steroid use, Prophylactic antibiotics use, Brittle diabetes, Ehler-Danlos syndrome, Osteoradionecrosis, Renal failure, organ transplantation, Anticoagulation therapy, unexplained hypersensitivity, Fibrous dysplasia, regional enteritis.
6. Diseases, or treatment that severely compromise healing, e.g., including radiation therapy.
7. Lack of adequate training of practitioner.
8. Poor patient motivation, such as Psychiatric disorders that interfere with patient understanding and compliance with necessary procedures, Unrealistic patient expectations, Unattainable prosthodontic reconstruction, Inability of patient to manage oral hygiene, Patient hypersensitivity to specific component of the procedure.
9. Electrosurgery: Dental implants are made of a metallic alloy; therefore, they are characterized by high conductivity. For this reason, electrosurgery is strictly contraindicated near dental implants.

#### Additional Possible Contraindications

Consider not using the device if the patient has one or more of the following medical situations: Chronic bleeding problems, psychological impairment, treatment with chemotherapeutic agents, metabolic bone or connective tissue diseases, treatment with corticosteroids, certain cardiac and vascular diseases, diabetes (uncontrolled), heavy smoker's usage, chronic renal disease, poor patient oral hygiene, bruxism, alcoholism.

#### Temporary Contraindications

1. Systemic infection, local oral and respiratory infection
2. Anatomical or Pathological implications such as.
  - a. Insufficient alveolar bone width and height to surround the implant with at least 1.5 millimeter of bone
  - b. Inadequate bone height where proper implant placement would

encroach within 2mm of the mandibular canal, sinus floor, etc.

3. Malignancies
4. Pregnancy
5. Inadequate bone volume unless augmentation procedure can be considered.

#### WARNINGS:

1. These devices are only to be used by trained professionals. The surgical and restorative techniques required to properly utilize these devices are highly specialized and complex procedures. Improper technique can lead to implant failure, loss of supporting bone, restoration fracture, screw loosening and aspiration. Therefore, no implantation should be performed without prior adequate training by a certified dentist.
2. With respect to pediatric patients, routine treatment is not recommended until the end of the jawbone growth has been properly documented.
3. This product is a single use device. Any attempt to reuse the device may cause severe risk to health such as product contamination, acute infection, degradation of mechanical properties leading to painful revision surgery for replacing the implant.
4. Do not use implants if the packaging is opened or damaged, as this can result with 004\_Instructions for Use CAD/CAM Abutments compromising sterility or device malfunction.

**Warning:** The AoN Implants Dental Implant/Abutment System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the AoN Implants Dental Implant/Abutment System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

1. **BREAKAGE** Implant and abutment fractures can occur when applied loads exceed the normal functional design tolerances of the implant components. Potential overloading conditions may result from deficiencies in implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 30 degrees, occlusal interferences causing excessive lateral forces, patient parafunction (e.g., bruxing clenching), improper casting procedures, inadequate prosthesis fit, and physical trauma.
2. Do not use excessive force exceeding insertion torque of 60Ncm: Overtightening an implant may lead to damage of the implant and/or internal connection, fracture or necrosis of bone site, compromising osseointegration.
3. Balanced occlusion.
4. Proper cementation.
5. Follow up of the patients.

#### CLEANING/STERILIZATION

The processes of cleaning and sterilization are intended to safeguard the health of patients and of all persons working in the surgery.

#### Cleaning and sterilization:

Cleaning: can be manually accomplished with hot water and a suitable neutral non-corrosive detergent, using plastic or nylon brushes (never steel wool or metal brushes) to remove any organic residues.

Subject AoN Implants			
Abutment	CYROTH Compatible	IS FOUR Compatible	KA-ROTH Compatible
Standard	Yes	Yes	Yes
Anatomic	Only Internal Hex	Yes	No
Multi-Unit	Yes	Yes	Yes
Ball Attachment	Yes	Yes	Yes
T-Base	Yes	Yes	Yes
T-Base CEREC	Yes	Yes	Yes
T-Base Multi-Unit	Yes	Yes	Yes
Conometric Prosthetic	Yes	Yes	Yes

The process of hand washing is effective if done with maximum care.

Cleaning can also be carried out in ultrasonic tanks, in a manner that is faster and safer, by fully immersing the pieces in washing liquid, and following the times recommended by the manufacturer of the equipment. Replace the washing solution according to the usage time limits recommended by the manufacturer. After ultrasonic cleaning, thoroughly rinse to remove the washing solution. Immediately after rinsing, before sterilization, inspect and dry them.


Sterilization: remember, first, that sterilization does not replace cleaning. Sterilization is a process that kills all microorganisms, and it must be undertaken with maximum care after cleaning. The sterilization methods most common and effective are the following:


1. Abutments, Prosthetics, CAD/CAM, Ball Attachments, Healing Caps, and Securement Screws - Autoclave in a Tuttnauer 387 (K003470) using Dynamic Air Removal (Pre-Vacuum) for 4 minutes at a temperature of 132°C with a Drying Time of 20 minutes.
2. Place the product before sterilization in a standard sterilization pouch (Steriking Wipak [K803293] or equivalent, compatible with steam sterilization process. Sets of instruments may be loaded into dedicated instrument trays or general-purpose trays for seam sterilization.
3. Only use FDA cleared sterilization packaging materials. Follow the autoclave manufacturer's instructions to sterilize the products. Care must be taken not to exceed the autoclave maximum recommended load. We recommend using a Tuttnauer 387 (or equivalent) and BIs K082756 or equivalent.

## 1.1 DISPOSAL

Dispose of according to facility standard operating procedures and protocols.

### EXPLANATION OF USED SYMBOLS

SYMBOL	STANDARD REFERENCE	SYMBOL TITLE	EXPLANATORY TEXT
	ISO 15223-1:2016 Reference no. 5.1.1. (ISO 7000-3082)	Manufacturer	Indicates the medical device manufacturer
	ISO 15223-1:2016 Reference no. 5.1.2	Authorized Representative in the European Community	Indicates the authorized representative in the European Community / European Union
	ISO 15223-1:2016 Reference no. 5.1.5. (ISO 7000-2492)	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified
	ISO 15223-1:2016 Reference no. 5.1.6. (ISO 7000-2493)	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified ISO 15223 Catalogue number ISO 7000 Catalog number
	ISO 15223-1:2016 Reference no. 5.2.7. (ISO 7000-2609)	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process

SYMBOL	STANDARD REFERENCE	SYMBOL TITLE	EXPLANATORY TEXT
	ISO 15223-1:2016 Reference no. 5.3.2. (ISO 7000-0624)	Keep away from sunlight	Indicates a medical device that needs protection from light sources
	ISO 15223-1:2016 and ISO 15223-1:2020(E) DRAFT Reference no. 5.3.4. (ISO 7000-0626)	Keep dry Keep away from rain	Indicates a medical device that needs protection from moisture ISO 15223 Keep dry ISO 7000 Keep away from rain
	ISO 15223-1:2016 and ISO 15223-1:2020 Reference no. 5.4.3. (ISO 7000-1641) Reference no. 5.4.3. (ISO 7000-1641)	Consult instructions for use Operator's manual; operating instructions	Indicates the need for the user to consult the instructions for use iso_15223 Consult instructions for use iso_grs_7000_1641 Operator's manual; operating instructions
	ISO 15223-1:2016 Reference no. 5.4.4. (ISO 7000-0434A)	Caution	To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action to avoid undesirable consequences
	ISO/DIS 15223-1:2020(E) DRAFT Reference no. 5.7.7	Medical device	Indicates the item is a medical device
	ASTM F2503 Reference no. Table 2, Symbol 7.3.3; 7.4.9.1; Fig. 9	(MR) Unsafe	3.1.14: An item which poses unacceptable risks to the patient, medical staff, or other persons within the MR environment
	EU 2017-745 EU 2017-746 Reference no. ANNEX V	CE marking	(43) 'CE marking of conformity' or 'CE marking' means a marking by which a manufacturer indicates that a device is in conformity with the applicable requirements set out in this Regulation and other applicable Union harmonization legislation providing for its affixing
	N/A	Prescription Use Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.
	N/A	Quantity	20.2 in the 2017/746 IVDR Regulation (j) where relevant, an indication of the net quantity of contents, expressed in terms of weight or volume, numerical count, or any combination of thereof, or other terms which accurately reflect the contents of the package.

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