## GENERAL ABUTMENTS

EN

### Instruction for Use

**AoN Implants US** prosthetic components: Z abutment, Multi Utility Temporary Abutment, Transepithelial Connector, Angle Overdenture Ball Abutment & OverDLock®

2 - 5

ES

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**AoN Implants US** componentes protésicos: Z Pilar, Multi Utility Pilar Temporal, Conector Transepithelial, Pilar de Bola Sobre Dentadura & OverDLock®

6 - 9

#### **Implant-Abutment Compatibility:**

The following identifies the intended combinations of abutment angulations and implant platforms to be used together:

- Conical Implants: Used with 15 and 25 degrees abutment angulations.
- Internal Hex Implants: Used with 10, 15, 20, 25 and 30 degrees abutment angulations.
- External Hex Implants: Used with 10 and 15 degrees abutment angulations.

#### **Warning:**

Small diameter implants and angled abutments are not recommended for use in the posterior region of the mouth.

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**Caution: Federal Law restricts this device to sale by or on the order of a licensed physician or health-care provider.**

## 1 Prosthetic components for dental implants

### 1.1 Description

The AoN implants' prosthetics called Healing Caps or Abutments in the variants: straight, angled, calcinable, temporary and customizable, are prosthetic structures, prefabricated for AoN Implant' dental implants.

### 1.2 Materials

All materials used are of high quality and conform to reference standards. AoN Implants' prosthetics are realized with materials supplied by qualified companies and conform to international reference standards. All materials used guarantee a high degree of biocompatibility. AoN Implants' prosthetics (including securement screws and healing caps) are sold in non-sterile packages. It is recommended, however, that the opening of packaging be done in aseptic surroundings and possibly over disinfected or sterile worktops and surfaces.

### 1.3 Indications for Use

Immediate loading when good primary stability is achieved and with appropriate occlusal loading for use in surgical (single-stage or two-stage procedures) and restorative applications for placement in the bone of the upper or lower jaw to provide support for prosthetic device such as artificial teeth and to restore the patient's chewing function. The System is intended to be used in either single tooth or multiple teeth applications. Short implants are to be used only with straight abutments.

### 1.4 Contraindications

The use of AoN Implants' prosthetics does not present any particular contraindications and it is suggested that you refer directly to the subsequent applications of dental implants and to their instructions to which it is necessary to refer for a clear indication.

Remember, however, that some contraindications are generally envisioned, such as: presence of chronic and acute infective diseases, systemic diseases, chronic sub-acute maxillary osteitis, diseases involving microvascular disorders.

The use of AoN Implants' prosthetics is discouraged in patients who do not cooperate, those who abuse of drugs and alcohol and those that have poor oral hygiene.

### 1.5 Warnings

The use of AoN Implants' prosthetics is reserved exclusively to personnel who are familiar with prosthetic techniques for dental implants.

The procedures to be employed are varied and complex; therefore, we defer to the competence of the specialist for the evaluation of the best technique to use.

The insertion of the prefabricated prosthetic structures on implants must be effectuated in surgically appropriate circumstances. An erroneous or improper use can cause the failure of the prosthetic operation.

## 1.6 Instructions for Use

Before insertion in the mouth, the Healing Caps, Securement Screws, and the abutments must be carefully cleaned and sterilized. On completion of use and when it is necessary to proceed with the disposal of the Healing Caps, Securement Screws, and abutments, in conventional ways, repeat the cleaning and sterilization.

The AoN Implants' prosthetics are of single use; this serves to reduce the risk of potential technical damage, due to previous uses, which could compromise the work or realizing the prosthetic operation.

**Healing cap:** on achieving osseointegration of the implant, in order to predispose the mucosa to a non-traumatically reception of the abutment, and to access the cover screw and replace it with the healing cap of a height adequate for the thickness of the mucosa, through a simple surgery. The healing time is around 15 to 30 days but may vary according to individual cases.

**Abutment:** upon healing of the mucosa, replace the Healing cap with a suitable abutment, fixing it to the bottom of the implant by means of the retention screw included in the pack, then secure the prosthesis on the abutment itself.

For the preparation of the prosthesis and the possible shape of the abutment, clearly it is necessary to take an imprint and execute precise work in the laboratory.

**Calcinable abutment:** complete the appropriate modeling working in the laboratory with the aid of the casts realized according to the imprints from the patient, and then proceed with melting and finishing of the abutment thus obtained that will be fixed to the implant by means of the retention screw which is provided.

Fix the prosthesis on the abutment. The retention screws of the abutment must be tightened with a force of less than 25N/m.

### 1.7 Use by Date

There are no special storage conditions although it is always preferable to store in a cool, dry place and not expose the product to direct sources of heat or sunlight.

The disposal of the product when no longer useable must follow the existing procedures on the subject of disposal of medical waste products.

## 1.8 Cleaning and Sterilization:

The processes of cleaning and sterilization are intended to safeguard the health of patients and of all persons working in the surgery.

**Cleaning:** can be manually accomplished with hot water and a suitable neutral non-corrosive detergent, using plastic or nylon brushes (never steel wool or metal brushes) to remove any organic residues.

The process of hand washing is effective if done with maximum care.

Cleaning can also be carried out in ultrasonic tanks, in a manner that is faster and safer, by fully immersing the pieces in washing liquid, and following the times

recommended by the manufacturer of the equipment. Replace the washing solution according to the usage time limits recommended by the manufacturer. After ultrasonic cleaning, thoroughly rinse to remove the

washing solution. Immediately after rinsing, before sterilization, inspect and dry them.

**Sterilization:** remember, first of all, that sterilization does not replace cleaning. Sterilization is a process that kills all microorganisms and it must be undertaken with maximum care after cleaning. The sterilization methods most common and effective are the following:

1. Abutments, Healing Caps, and Securement Screws - Autoclave in a Tuttnauer 387 (K003470) using Dynamic Air Removal (Pre-Vacuum) for 20 minutes at a temperature of 132°C with a Drying Time of 20 minutes.
2. Place the product before sterilization in a standard sterilization pouch (Striking Wipak [K803293] or equivalent, compatible with steam sterilization process. Sets of instruments may be loaded into dedicated instrument trays or general-purpose trays for steam sterilization.
3. Only use FDA cleared sterilization packaging materials. Follow the autoclave manufacturer's instructions to sterilize the products. Care must be taken not to exceed the autoclave maximum recommended load. We recommend using a Tuttnauer 387 (or equivalent) and BIs K082756 or equivalent.

## 1.9 Instructions for Use Angled Overdenture Ball Abutment (Sphero Block)

### DESCRIPTION

The Angled Overdenture Ball Abutment is a universal hinge, resilient attachment for endosseous implants and angled or straight multi-unit abutments in the mandible or maxilla that is designed to restore masticatory function. The attachment system is patient-removable and therefore allows for the prosthesis to be removed and replaced by the patient.

The Angled Overdenture Ball Abutment is designed for use with full arch overdentures or partial dentures retained in whole or in part by endosseous implants in the mandible or maxilla.

The Angled Overdenture Ball Abutment is not recommended for use with a single implant with divergence of greater than 10 degrees from vertical. Not appropriate where a totally rigid connection is required.

Angled Overdenture Ball Abutment (Sphero Block)	
Description	The Angled Overdenture Ball Abutment is an abutment with an internal and external hex that connects to the reference K192015 implants for use in edentulous situations
Restoration Technique	to be screwed on the implant
Characteristics	The abutment is intended for single use. It is fixated using a hex driver that turns and tightens the securement screws to the implant. The prosthetic is not used with another prosthetic.
Family Type	Overdenture Prosthetic
Implant Attachment	Securement Screw
Prosthetic Attachment	Retention with Cap

## 1.10 Instructions for Use. OverDLock Attachment System

The OverDLock Attachment System includes: Dental Implant Abutments, Denture Attachment Housing, Retention Inserts, Ancillary

Processing Components (including: abutment analogs, processing spacer, impression coping, block-out spacer and processing insert), and Tooling (including: Retention Insert Tool and abutment drivers for torque wrenches and latch handpieces).

This document contains the most current Instructions for Use. Please read and retain.

The OverDLock Attachment System is a universal hinge, resilient attachment for endosseous implants and angled or straight multi-unit abutments in the mandible or maxilla that is designed to restore masticatory function. The attachment system is patient-removable and therefore allows for the prosthesis to be removed and replaced by the patient.

The OverDLock Attachment System is designed to facilitate patient removal of a dental prosthesis for use with full arch overdentures or partial dentures retained in whole or in part by endosseous implants in the mandible or maxilla.

Not recommended for use with a single implant with divergence of greater than 20 degrees from vertical. Not appropriate where a totally rigid connection is required.

OverDLock Attachment Implant Compatibility	
Implant	Configurations
AoN Implants Dental Implants Cyroth and IS-Four (K192015)	3.5mm Platform (internal hex), 2.8mm (conical connection) and 3.4/4.1mm (external hex)

### 1.10.1 Warnings and Precautions

Product should be inspected for integrity prior to use. Product from damaged packaging should not be used on patients. In the event that the packaging is damaged, the damaged packaging with the product should be returned to the manufacturer and a replacement will be provided only if damage to packaging is caused by product shipment.

As surgical instruments are susceptible to damage and wear, they should be inspected before each use. Markings should be visible and legible. Any reusable instrument should be replaced if damage or wear is present to ensure proper functionality. The number of uses will vary and depends on a variety of factors including but not limited to bone density encountered, handling, proper cleaning, autoclave exposure, and storage conditions (do not store tools or instruments wet). Over time, repeat sterilization may affect appearance and visibility of markings. When applicable to the surgical instrument, check the connection feature for wear to ensure the connection is not damaged.

Patient evaluation including the determination of the general health, oral hygiene habits and status, motivation towards good dental care, and anatomical acceptability prior to the placement of the implant attachments as part of restorative process, is critical. Thorough evaluation of the patient's medical status and health history is mandatory. Treatment planning is vital to the success of the implant and prosthesis.

The use of this attachment system requires that the clinician be thoroughly familiar with the product and the method for its use and application. The clinician must also use

reasonable judgment in deciding when and where to use the product.

The OverDLock Attachment System components are single-use devices and are provided non-sterile. Single-use devices must not be re-used or re-sterilized. Reuse of a single-use device may cause harm to the patient in the transfer of blood, tissue or saliva that may contain infectious disease. Single-use devices may not function as intended if re-sterilized and may result in an improper surgical procedure and lead to improper function or failure of the device.

OverDLock Retention Inserts: The inadvertent re-use of OverDLock nylon Retention Inserts could cause loss of retention of the overdenture due to wear from previous use or damage during removal.

The surgical instruments and tools of the OverDLock Attachment System are multi-use devices. Reusable tools and instruments must be cleaned and sterilized prior to reuse on patients.

Tooling: The OverDLock Tools are designed for multiple uses and are provided NON-STERILE. Follow the instructions provided here within for proper sterilization of non-sterile components and the instructions for cleaning and re-sterilization processes of reusable components.

The OverDLock Attachment System restorative components and instruments are supplied NON-STERILE and should be sterilized prior to use on patients.

### 1.10.2 Sterilization

Inserts, assemblies which include inserts (such as the denture housing assembly and impression coping), and the blockout spacer are heat sensitive. An FDA approved liquid chemical sterilant is recommended for use with these components to achieve appropriate disinfection/sterilization of these components. These components may be sterilized/disinfected using a liquid chemical sterilant according to the manufacturer's instructions.

Abutments and instruments, including drivers and the system tool must be sterilized prior to use on patients. These devices may be sterilized by autoclave sterilization using the following parameters. For gravity cycle, place components in autoclave bag; and for Pre-Vacuum Cycle, wrap the component with autoclave wrap material and secure wrap with autoclave tape. Wrap the components using a wrap that is FDA-cleared for the indicated cycles.

Autoclave Sterilization Parameters are listed below:

Part Description	Cycle Type	Temperature	Exposure Time	Drying Time
OverDLock Abutment and Instruments including System Tool and Abutment Drivers	Gravity	132°C / 270°F	15 Minutes	30 Minutes
	Pre-Vacuum	132°C / 270°F	4 Minutes	20 Minutes

Re-sterilizable instruments should be dried completely and stored in a clean and dry location at normal room temperature. Prior to instrument use, the exterior of any sterilized packaging should be inspected for integrity. Care must be exercised in the handling of wrapped or autoclave bagged instrument kits or instruments to prevent damage to the sterile barrier. If damage to the sterile barrier is

observed, re-sterilization is recommended for reusable devices only. Single Use devices should not be re-sterilized.

### 1.10.3 Cleaning

Reusable instruments must be cleaned in between patient uses according to the following instructions.

- Disconnect or disassemble the instruments.
- Soak the instruments in enzymatic cleaning solution mixed according to manufacturer's instructions by completely submerging them for 20 minutes. Scrub the components using a soft-bristled, nylon brush until all soil has been removed.
- Remove the instruments from the enzymatic cleaning solution and rinse in tap water for a minimum of 3 minutes. Make sure to thoroughly flush internal holes/crevices of the instruments that have difficult to reach areas.
- Visually inspect instruments and tools for cleanliness and presence of residual debris. If additional cleaning is needed, place instruments in ultrasonic cleaner with enzymatic cleaning solution prepared according to manufacturer's instructions making sure that they are completely submerged and sonicate for 10 minutes.
- Remove the instruments from the ultrasonic cleaner, and rinse for 3 minutes making sure to thoroughly flush cleaning solution out of the holes/crevices and/or difficult to reach areas.
- Remove excess moisture from the instruments with a clean, absorbent, non-shedding wipe and allow instruments to dry completely.

### 1.10.4 Healing Phase

For delayed loading protocols: Relieve the denture to ensure the abutments are not in contact with any denture acrylic. A soft liner may be added to the denture to ensure patient comfort during the healing phase.

### 1.10.5 Patient Care

Good oral hygiene is vital to success with the OverDLock Attachment System. The patient should be made aware of the following:

- Attachments must be thoroughly cleaned each day to prevent buildup of plaque biofilm. The patient should use a soft nylon bristle or end-tufted toothbrush with a non-abrasive toothpaste to clean the Abutments.
- The coarse particles in abrasive toothpastes may scratch the surfaces of the attachments and cause additional plaque accumulation.
- An irrigation system is recommended to flush out debris from the inside of the OverDLock Retention Inserts.
- OverDLock Retention Inserts are made of a soft plastic material (nylon) to allow the overdentures to be removed and replaced regularly. Plastic materials are subject to wear as part of normal use and may require replacement.
- Bruxism may reduce the longevity of Retention Inserts.

Patients should be instructed to maintain routine follow-up visits for hygiene and evaluation of attachment function. Should a patient experience any discomfort or loss of

overdenture retention, they should consult a dental professional. Follow-up visits are recommended at 6-month intervals.

#### Inserting and Removing Overdentures

The patient should be instructed on how to properly insert the Overdenture. The patient should make sure they can feel that it is positioned over the Abutments prior to applying pressure. The patient should use both hands and press down on each side and firmly snap the Overdenture into place.

**NOTE:** The patient must NOT bite their Overdentures into place, as this force will result in improper wear of the Abutments and Retention Inserts. The Overdenture is to be removed by the patient by placing their thumbs under the edges of the Overdenture flanges and pulling each side upward (mandibular denture) or downward (maxillary denture) simultaneously. Use of the tongue may aid in removal. Once removed, thorough cleaning is recommended.

#### 1.10.6 Cleaning of Implant-Retained Overdentures

Instruct the patient to follow the protocol below to ensure the longevity of their Overdenture.

1. Fill a washing basin with warm water to prevent fracture of the Overdenture. Apply non-abrasive toothpaste onto a soft nylon bristle or end-tufted toothbrush and thoroughly clean every surface of the Overdenture.
2. Each night removes the Overdenture and rinse with plain water.

#### 1.10.7 Additional Information

Traditional restorative protocols should be followed to process the attachments into the patient's Overdenture. Standard Overdenture care and maintenance should be followed in order to ensure the longevity of each restoration.

<b>OverDLock</b>	
Description	The OverDLock is an abutment with an internal hex, conical, or external hex platform that connects to the reference K192015 implants for use in edentulous situations
Restoration Technique	to be screwed on the implant
Characteristics	The OverDLock is designed to retain overdentures or partial dentures in the mandible or maxilla. The abutment is intended for single use. It is fixated using a hex driver that turns and tightens directly to the implant. The prosthetic is not used with another prosthetic.
Family Type	Overdenture Prosthetic
Implant Attachment	Direct Screwed
Prosthetic Attachment	Retention with Cap

<b>OverDLock Cap Assembly</b>	
Description	The OverDLock Cap is intended to be used in conjunction with OverDLock Retention Insert the and OverDLock abutment that connects to the reference K192015 implants for use in edentulous situations.
Restoration Technique	to be pressed on the OverDLock abutment together with OverDLock Retention Insert

<b>OverDLock Cap Assembly</b>	
Characteristics	The OverDLock Cap is designed to retain overdentures or partial dentures in the mandible or maxilla. The cap is intended for single use. The prosthetic is not used with another prosthetic.
Family Type	Overdenture Prosthetic
Prosthetic Attachment	Retention with Cap
Sterilization	No - Sterile

<b>OverDLock Retention Inserts</b>	
Description	The OverDLock Retention Inserts are intended to be used in conjunction with the OverDLock Caps and OverDLock abutments that connects to the reference K192015 implants for use in edentulous situations. The component is available in various hardness's according to the prosthetist's requests. Hardness varies according to the color of the plastic component.
Restoration Technique	to be attached between the OverDLock abutment and OverDLock cap.
Characteristics	The OverDLock Retention Insert is designed to retain overdentures or partial dentures in the mandible or maxilla. The insert is intended for single use. The prosthetic is not used with another prosthetic.
Family Type	Overdenture Prosthetic
Prosthetic Attachment	Retention with Cap
Material	Nylon
Colors	Blue, Pink, White, Purple, Gray, Red, Yellow, Green
Sterilization	No - Sterile

<b>OverDLock Processing Insert</b>	
Description	The OverDLock abutment is coupled with this plastic component which is fixed in the prosthesis by means of a titanium cap. The OverDLock Processing Insert is used for the processing of the denture cap into the denture.
Insert Pivot Types	Standard, Standard PVD Cap, Extended, Standard Cap
Insert Retention Type	NA
Restoration Technique	to be attached between the OverDLock abutment and OverDLock cap.
Characteristics	It is used for the process of gluing the prosthesis on the abutment in the patient's mouth. It does not generate retention in the abutment. It is used in the mouth for no more than 5-10 minutes at a time.
Family Type	OverDLock Attachment System
Prosthetic Attachment	Retention with Cap
Material	Low Density Polyethylene
Colors	Black
Sterilization	No - Sterile

<b>OverDLock Processing Spacer</b>	
Description	The OverDLock abutment is coupled with this plastic component which is fixed in the prosthesis by means of a titanium cap. The OverDLock Processing Spacer is not used in the mouth at all – it is a part that is designed for a laboratory technician to make the proper space in the bottom of a denture to accommodate the denture cap.
Restoration Technique	NA
Characteristics	Laboratory piece
Family Type	OverDLock Attachment System
Prosthetic Attachment	Retention with Cap
Material	Delrin
Colors	White
Sterilization	No - Sterile



<b>OverDLock Block-Out Spacer</b>	
Description	The OverDLock abutment is coupled with this plastic component which is fixed in the prosthesis by means of a titanium cap. The OverDLock Block-Out Spacer is only used in the mouth during the procedure to pick up the denture caps into the bottom of the denture.
Restoration Technique	to be attached between the OverDLock abutment and OverDLock cap.
Characteristics	It is used in the mouth for no more than 5-10 minutes at a time.
Family Type	OverDLock Attachment System
Prosthetic Attachment	Retention with Cap
Material	Silicone
Colors	White
Sterilization	No - Sterile

#### 1.11 DISPOSAL

Dispose of according to facility standard operating procedures and protocols.

#### EXPLANATION OF USED SYMBOLS

SYMBOL	STANDARD REFERENCE	SYMBOL TITLE	EXPLANATORY TEXT
	ISO 15223-1:2016 Reference no. 5.1.1. (ISO 7000-3082)	Manufacturer	Indicates the medical device manufacturer
	ISO 15223-1:2016 Reference no. 5.1.2	Authorized Representative in the European Community	Indicates the authorized representative in the European Community / European Union
	ISO 15223-1:2016 Reference no. 5.1.5. (ISO 7000-2492)	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified
	ISO 15223-1:2016 Reference no. 5.1.6. (ISO 7000-2493)	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified ISO 15223 Catalogue number ISO 7000 Catalog number
	ISO 15223-1:2016 Reference no. 5.2.7. (ISO 7000-2609)	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process
	ISO 15223-1:2016 Reference no. 5.3.2. (ISO 7000-0624)	Keep away from sunlight	Indicates a medical device that needs protection from light sources
	ISO 15223-1:2016 and ISO 15223-1:2020(E) DRAFT Reference no. 5.3.4. (ISO 7000-0626)	Keep dry Keep away from rain	Indicates a medical device that needs protection from moisture ISO 15223 Keep dry ISO 7000 Keep away from rain
	ISO 15223-1:2016 and ISO 15223-1:2020 Reference no. 5.4.3. (ISO 7000-1641) Reference no. 5.4.3. (ISO 7000-1641)	Consult instructions for use Operator's manual; operating instructions	Indicates the need for the user to consult the instructions for use iso_15223 Consult instructions for use iso_grs_7000_1641 Operator's manual; operating instructions

SYMBOL	STANDARD REFERENCE	SYMBOL TITLE	EXPLANATORY TEXT
	ISO 15223- 1:2016 Reference no. 5.4.4. (ISO 7000-0434A)	Caution	To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action to avoid undesirable consequences
	ISO/DIS 15223-1:2020(E) DRAFT Reference no. 5.7.7	Medical device	Indicates the item is a medical device
	ASTM F2503 Reference no. Table 2, Symbol 7.3.3; 7.4.9.1; Fig. 9	(MR) Unsafe	3.1.14: An item which poses unacceptable risks to the patient, medical staff, or other persons within the MR environment
	EU 2017-745 EU 2017-746 Reference no. ANNEX V	CE marking	(43) 'CE marking of conformity' or 'CE marking' means a marking by which a manufacturer indicates that a device is in conformity with the applicable requirements set out in this Regulation and other applicable Union harmonization legislation providing for its affixing
	N/A	Prescription Use Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.
	N/A	Quantity	20.2 in the 2017/746 IVDR Regulation (j) where relevant, an indication of the net quantity of contents, expressed in terms of weight or volume, numerical count, or any combination of thereof, or other terms which accurately reflect the contents of the package.

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